

prodoctor®

Erken Gebelik Testi | Early Pregnancy Test

EARLY PREGNANCY TEST

99.8% ACCURATE
RESULT



EARLY
PRACTICAL
FAST
HYGIENIC

PROFESSIONAL SERIES

A 'missed period' refers to the day after the expected period
75% of positive (pregnant) results can be detected 7 days before





ABOUT



PRODOCTOR

Prodoctor is the best-selling home pregnancy and fertility test in Turkey. Consumers prefer the Prodoctor brand because it delivers the accurate information they trust. We owe this trust to our professional and experienced teammates . Prodoctor continues to give consumers the fastest, practical and accurate results with its experience of more than 15 years and 5 different product types. To contact Prodoctor's professional sales team.



HOW PRODOCTOR DETECS PREGNANCY?

Over **99%** accurate at detecting pregnancy from the day of the expected period.

Prodoctor early pregnancy test detects the level of human chorionic gonadotrophin (hCG) hormone in the urinary. This hormone produced by the embryo during the earliest stages of pregnancy.

Levels of hCG in serum and urine rise rapidly during the first days of pregnancy by approximately 50% per day. In the first 3 weeks of pregnancy the hCG hormone level can be measured at different values, but the hCG hormone levels after 4 weeks are quite similar. Using the Prodoctor after the first week is the ideal times for accurate and early pregnancy detection for women.

If used on or after the day of the expected period, the Prodoctor Early Pregnancy Test is more than 99% accurate at detecting pregnancy.

	1-7 Days	8-14 Days	15-21 Days	22-28 Days	28 + Days
0 mIU/mL	+				
25mIU/mL	+				
100mIU/mL	+	+			
300mIU/mL		+	+		
1000mIU/mL			+	+	
10000mIU/mL +				+	+

Median hCG concentration (mIU/mL)

Duration of pregnancy (days)

TEST

● **PRODOCTOR EARLY PREGNANCY TEST**


Prodoctor Early Pregnancy Test Features With Ergonomic Design

- A long and lightweight design crafted for comfortable and easy use
- A wide cottony structure designed to obtain an adequate sample from the urinary
- Double layer gelatin coating prepared for hygienic and reliability
- A large easy-to-read result window

Prodoctor early pregnancy test helps women to use correctly and practically.



**PROVIDES
EASY-TO-READ
RESULT**



CONTROL RESULT

Two lines indicates a positive (pregnant) result. The appearance of the line in the result window may vary; however, it is still a 'positive' result no matter how faint or dark the line is.



CONTROL RESULT

A line only in the control indicates a negative result at the time of testing.

- The test screen shows results between 3 and 10 minutes. If no line has appeared in the control window after 10 minutes, the test has not worked.
- If the user is testing before her expected period, the level of hCG may not yet be high enough to be detected and she is advised to test again when her period is expected.
- Ectopic pregnancy and ovarian cysts can give misleading results.


PLUS

● PRODOCTOR PLUS EARLY PREGNANCY TEST


Prodoctor Plus with weekly pregnancy test.

- Wide and large window designed easy-to-read
- Improved tip for accurate and superior fast reading
- More precise reading with its innovative and technology for weekly pregnancy detection


The prodoctor plus early pregnancy test helps women determine what week of pregnancy they are in.




NEGATIVE




1 WEEK



2 WEEK



3 WEEK



4 + WEEK

NOTE

- In order for the test to give more accurate results, it is more accurate to test with the first urine of the day. Consuming fluids before using the test will dilute the hCG hormone in the urine, making it difficult to determine the correct week.
- After using the test, if one or both of the 2 control lines do not appeared, the test is considered invalid.
- Testing within 14 days of using fertility meds and birth control pills can affect the result.

PRODOCTOR EARLY PREGNANCY TESTS ••

EARLY

Fastest and Earliest Pregnancy Test

SENSITIVE

Can be used 7 days before the missed period

ACCURATE

Gives women the most accurate results with 99.8% accuracy

RELIABLE

Turkey's most preferred brand

PRACTICAL

Easy and ergonomic design

prodoctor®

Erken Gebelik Testi | Early Pregnancy Test

www.prodoctor.com.tr

+90 532
577 75 50

Işıklar Mahallesi Hilmedede Sokak
2/B SELÇUKLU / KONYA

info
@prodoctor.com.tr